

Product Datasheet

Human GZMK ELISA Kit (orb1200563)

Description	Human Granzyme K (GZMK) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49863
Sensitivity	13.1pg/mL
Entrez	3003
Expiration Date	Please enquire.