

## Product Datasheet

### Human GLUT1 ELISA Kit (orb1200586)

|                            |                                                                 |
|----------------------------|-----------------------------------------------------------------|
| <b>Description</b>         | Human Glucose Transporter 1 (GLUT1)<br>ELISA Kit                |
| <b>Reactivity</b>          | Human                                                           |
| <b>Range</b>               | 0.312-20ng/mL                                                   |
| <b>Tested Applications</b> | ELISA                                                           |
| <b>Note</b>                | For research use only                                           |
| <b>Application notes</b>   | sandwich                                                        |
| <b>Sample Types</b>        | tissue homogenates, cell lysates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <b>P11166</b>                                                   |
| <b>Sensitivity</b>         | 0.123ng/mL                                                      |
| <b>Entrez</b>              | <b>6513</b>                                                     |
| <b>Expiration Date</b>     | Please enquire.                                                 |