

Product Datasheet

Human CYSLTR2 ELISA Kit (orb1200655)

Description	Human Cysteinyl Leukotriene Receptor 2 (CYSLTR2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9NS75
Sensitivity	0.108ng/mL
Entrez	57105
Expiration Date	Please enquire.