

Product Datasheet

Human CX3CR1 ELISA Kit (orb1200660)

Description	Human Chemokine C-X3-C-Motif Receptor 1 (CX3CR1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P49238
Sensitivity	0.055ng/mL
Entrez	1524
Expiration Date	Please enquire.