

Product Datasheet

Human APOA5 ELISA Kit (orb1200757)

Description	Human Apolipoprotein A5 (APOA5) ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6Q788
Sensitivity	26.73pg/mL
Entrez	116519
Expiration Date	Please enquire.