

Product Datasheet

Human ALOX5 ELISA Kit (orb1200772)

Description	Human Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09917
Sensitivity	0.057ng/mL
Entrez	240
Expiration Date	Please enquire.