

Product Datasheet

Human ALDH7A1 ELISA Kit (orb1200775)

Description	Human Aldehyde Dehydrogenase 7 Family, Member A1 (ALDH7A1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49419
Sensitivity	0.29ng/mL
Entrez	501
Expiration Date	Please enquire.