

Product Datasheet

Human ADRP ELISA Kit (orb1200782)

Description	Human Adipose Differentiation Related Protein (ADRP) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99541
Sensitivity	0.057ng/mL
Entrez	123
Expiration Date	Please enquire.