

## Product Datasheet

### Human ACTb ELISA Kit (orb1200795)

|                            |                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------|
| <b>Description</b>         | Human Actin Beta (ACTb) ELISA Kit                                                     |
| <b>Reactivity</b>          | Human                                                                                 |
| <b>Range</b>               | 1.56-100ng/mL                                                                         |
| <b>Tested Applications</b> | ELISA                                                                                 |
| <b>Note</b>                | For research use only                                                                 |
| <b>Application notes</b>   | sandwich                                                                              |
| <b>Sample Types</b>        | tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>P60709</b>                                                                         |
| <b>Sensitivity</b>         | 0.63ng/mL                                                                             |
| <b>Entrez</b>              | <b>60</b>                                                                             |
| <b>Expiration Date</b>     | Please enquire.                                                                       |