

Product Datasheet

Rat IL20 ELISA Kit (orb1200816)

Description	Rat Interleukin 20 (IL20) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	D4A2T8
Sensitivity	5.7pg/mL
Entrez	681112
Expiration Date	Please enquire.