

Product Datasheet

Bovine PLGF ELISA Kit (orb1200840)

Description	Bovine Placenta Growth Factor (PLGF) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XS47
Sensitivity	6.8pg/mL
Entrez	280894
Expiration Date	Please enquire.