

Product Datasheet

Bovine IL2 ELISA Kit (orb1200852)

Description	Bovine Interleukin 2 (IL2) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05016
Sensitivity	5.9pg/mL
Entrez	280822
Expiration Date	Please enquire.