

Product Datasheet

Bovine IL18 ELISA Kit (orb1200859)

Description	Bovine Interleukin 18 (IL18) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9TU73
Sensitivity	6.1pg/mL
Entrez	281249
Expiration Date	Please enquire.