

Product Datasheet

Bovine IL17 ELISA Kit (orb1200861)

Description	Bovine Interleukin 17 (IL17) ELISA Kit
Reactivity	Bovine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	B3GDZ3
Sensitivity	11.2pg/mL
Entrez	282863
Expiration Date	Please enquire.