

Product Datasheet

Porcine IL16 ELISA Kit (orb1200862)

Description	Porcine Interleukin 16 (IL16) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q76KI7
Sensitivity	8.68pg/mL
Entrez	396589
Expiration Date	Please enquire.