

Product Datasheet

Porcine IL15 ELISA Kit (orb1200864)

Description	Porcine Interleukin 15 (IL15) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q95253
Sensitivity	12.8pg/mL
Entrez	397683
Expiration Date	Please enquire.