

Product Datasheet

Bovine IL12A ELISA Kit (orb1200868)

Description	Bovine Interleukin 12A (IL12A) ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Uniprot ID	P54349
Sensitivity	5.9pg/mL
Entrez	281856
Expiration Date	Please enquire.