

Product Datasheet

Bovine IGFBP4 ELISA Kit (orb1200874)

Description	Bovine Insulin Like Growth Factor Binding Protein 4 (IGFBP4) ELISA Kit
Reactivity	Bovine
Range	93.7-6000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q05716
Sensitivity	45.87pg/mL
Entrez	282262
Expiration Date	Please enquire.