

Product Datasheet

Porcine IFN γ ELISA Kit (orb1200885)

Description	Porcine Interferon Gamma (IFN γ) ELISA Kit
Reactivity	Porcine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17803
Sensitivity	5.9pg/mL
Entrez	396991
Expiration Date	Please enquire.