

Product Datasheet

Bovine FGF2 ELISA Kit (orb1200896)

Description	Bovine Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	F1N3T9
Sensitivity	6.5pg/mL
Entrez	281161
Expiration Date	Please enquire.