

Product Datasheet

Bovine CTGF ELISA Kit (orb1200904)

Description	Bovine Connective Tissue Growth Factor (CTGF) ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O18739
Sensitivity	30.23pg/mL
Entrez	281103
Expiration Date	Please enquire.