

Product Datasheet

Rat TFF2 ELISA Kit (orb1200959)

Description	Rat Trefoil Factor 2 (TFF2) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q09030
Sensitivity	0.064ng/mL
Entrez	116592
Expiration Date	Please enquire.