

Product Datasheet

Mouse SPA ELISA Kit (orb1200968)

Description	Mouse Surfactant Associated Protein A (SPA) ELISA Kit
Reactivity	Mouse
Range	3.12-200pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates or other biological fluids.
Uniprot ID	P35242
Sensitivity	1.45pg/mL
Entrez	20387
Expiration Date	Please enquire.