

Product Datasheet

Mouse SOD3 ELISA Kit (orb1200969)

Description	Mouse Superoxide Dismutase 3, Extracellular (SOD3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	009164
Sensitivity	0.064ng/mL
Entrez	20657
Expiration Date	Please enquire.