

Product Datasheet

Mouse PON1 ELISA Kit (orb1200997)

Description	Mouse Paraoxonase 1 (PON1) ELISA Kit
Reactivity	Mouse
Range	1.25-80ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P52430
Sensitivity	0.51ng/mL
Entrez	18979
Expiration Date	Please enquire.