

Product Datasheet

Mouse PKR1 ELISA Kit (orb1200998)

Description	Mouse Prokineticin Receptor 1 (PKR1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9JKL1
Sensitivity	0.04ng/mL
Entrez	58182
Expiration Date	Please enquire.