

Product Datasheet

Mouse ON ELISA Kit (orb1201027)

Description	Mouse Osteonectin (ON) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07214
Sensitivity	0.144ng/mL
Entrez	20692
Expiration Date	Please enquire.