

Product Datasheet

Mouse NSE ELISA Kit (orb1201035)

Description	Mouse Enolase, Neuron Specific (NSE) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17183
Sensitivity	0.069ng/mL
Entrez	13807
Expiration Date	Please enquire.