

Product Datasheet

Mouse MYH1 ELISA Kit (orb1201050)

Description	Mouse Myosin Heavy Chain 1, Skeletal Muscle, Adult (MYH1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q5SX40
Sensitivity	7.3pg/mL
Entrez	17879
Expiration Date	Please enquire.