

Product Datasheet

Mouse MMP7 ELISA Kit (orb1201057)

Description	Mouse Matrix Metalloproteinase 7 (MMP7) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q10738
Sensitivity	26.5pg/mL
Entrez	17393
Expiration Date	Please enquire.