

Product Datasheet

Mouse MMP2 ELISA Kit (orb1201059)

Description	Mouse Matrix Metalloproteinase 2 (MMP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P33434
Sensitivity	0.065ng/mL
Entrez	17390
Expiration Date	Please enquire.