

Product Datasheet

Mouse MK ELISA Kit (orb1201065)

Description	Mouse Midkine (MK) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12025
Sensitivity	0.059ng/mL
Entrez	17242
Expiration Date	Please enquire.