

Product Datasheet

Mouse MIP3b ELISA Kit (orb1201066)

Description	Mouse Macrophage Inflammatory Protein 3 Beta (MIP3b) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	070460
Sensitivity	5.9pg/mL
Entrez	24047
Expiration Date	Please enquire.