

Product Datasheet

Mouse LIPD ELISA Kit (orb1201080)

Description	Mouse Lipase, Lipoprotein (LIPD) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11152
Sensitivity	24pg/mL
Entrez	16956
Expiration Date	Please enquire.