

Product Datasheet

Mouse LAT ELISA Kit (orb1201086)

Description	Mouse Linker For Activation Of T-Cell (LAT) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O54957
Sensitivity	0.039ng/mL
Entrez	16797
Expiration Date	Please enquire.