

Product Datasheet

Mouse ICAM5 ELISA Kit (orb1201102)

Description	Mouse Intercellular Adhesion Molecule 5 (ICAM5) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q60625
Sensitivity	0.57ng/mL
Entrez	15898
Expiration Date	Please enquire.