

Product Datasheet

Mouse GFAP ELISA Kit (orb1201137)

Description	Mouse Glial Fibrillary Acidic Protein (GFAP) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P03995
Sensitivity	10.0pg/mL
Entrez	14580
Expiration Date	Please enquire.