

Product Datasheet

Mouse GAL8 ELISA Kit (orb1201143)

Description	Mouse Galectin 8 (GAL8) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9JL15
Sensitivity	0.055ng/mL
Entrez	56048
Expiration Date	Please enquire.