

Product Datasheet

Mouse GAL7 ELISA Kit (orb1201144)

Description	Mouse Galectin 7 (GAL7) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9CRB1
Sensitivity	0.049ng/mL
Entrez	16858
Expiration Date	Please enquire.