

Product Datasheet

Mouse FGL2 ELISA Kit (orb1201156)

Description	Mouse Fibrinogen Like Protein 2 (FGL2) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P12804
Sensitivity	5.4pg/mL
Entrez	14190
Expiration Date	Please enquire.