

Product Datasheet

Mouse CP ELISA Kit (orb1201194)

Description	Mouse Ceruloplasmin (CP) ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q61147
Sensitivity	0.27ng/mL
Entrez	12870
Expiration Date	Please enquire.