

Product Datasheet

Mouse CDHE ELISA Kit (orb1201204)

Description	Mouse Cadherin, Epithelial (CDHE) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P09803
Sensitivity	0.078ng/mL
Entrez	12550
Expiration Date	Please enquire.