

Product Datasheet

Mouse CD30L ELISA Kit (orb1201206)

Description	Mouse Cluster Of Differentiation 30 Ligand (CD30L) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P32972
Sensitivity	0.032ng/mL
Entrez	21949
Expiration Date	Please enquire.