

Product Datasheet

Mouse CCND1 ELISA Kit (orb1201208)

Description	Mouse Cyclin D1 (CCND1) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P25322
Sensitivity	5.7pg/mL
Entrez	12443
Expiration Date	Please enquire.