

Product Datasheet

Mouse ATP1b3 ELISA Kit (orb1201231)

Description	Mouse ATPase, Na ⁺ /K ⁺ Transporting Beta 3 Polypeptide (ATP1b3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P97370
Sensitivity	0.066ng/mL
Entrez	11933
Expiration Date	Please enquire.