

Product Datasheet

Mouse AQP4 ELISA Kit (orb1201234)

Description	Mouse Aquaporin 4 (AQP4) ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P55088
Sensitivity	23.7pg/mL
Entrez	11829
Expiration Date	Please enquire.