

Product Datasheet

Mouse aHSG ELISA Kit (orb1201252)

Description	Mouse Alpha-2-Heremans Schmid Glycoprotein (aHSG) ELISA Kit
Reactivity	Mouse
Range	9.37-600ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P29699
Sensitivity	3.64ng/mL
Entrez	11625
Expiration Date	Please enquire.