

## Product Datasheet

### Mouse ABCG2 ELISA Kit (orb1201260)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | Mouse ATP Binding Cassette Transporter G2 (ABCG2) ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                       |
| <b>Range</b>               | 0.156-10ng/mL                                               |
| <b>Tested Applications</b> | ELISA                                                       |
| <b>Note</b>                | For research use only                                       |
| <b>Application notes</b>   | sandwich                                                    |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.              |
| <b>Uniprot ID</b>          | <b>Q7TMS5</b>                                               |
| <b>Sensitivity</b>         | 0.055ng/mL                                                  |
| <b>Entrez</b>              | <b>26357</b>                                                |
| <b>Expiration Date</b>     | Please enquire.                                             |