

Product Datasheet

Mouse a1AGP ELISA Kit (orb1201265)

Description	Mouse Alpha-1-Acid Glycoprotein (a1AGP) ELISA Kit
Reactivity	Mouse
Range	62.5-4000ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	J3JRU4
Sensitivity	19.3ng/mL
Entrez	18407
Expiration Date	Please enquire.