

Product Datasheet

Human TLR9 ELISA Kit (orb1201309)

Description	Human Toll Like Receptor 9 (TLR9) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NR96
Sensitivity	0.059ng/mL
Entrez	54106
Expiration Date	Please enquire.